

Wiping Cloth Sanitizer Buckets



- Store wet wiping cloths in a sanitizer bucket between uses.
- Wiping cloths used with raw animal products must be kept separate from cloths used for other purposes.
- Check sanitizer strength with the proper sanitizer test kit:
 - Chlorine: 100 ppm
 - Quaternary Ammonium: 200 ppm
- Change sanitizing solution every 2 hours or as needed.



1-3 West Winter St. ♦ P.O. Box 570 ♦ Delaware, Ohio 43015
www.delawarehealth.org ♦ Fax 740-368-1736 ♦ Phone 740-368-1700

The Delaware General Health District operates in accordance with Title VI of the Civil Rights Act of 1964.